

Qdoser^{GO}
Qdoser^{ONE}



AUTOMATIC DOSING UNIT FOR MANUAL GRINDING & POLISHING MACHINES

QDOSER GO

The Qdoser GO automatic dosing unit may be used in combination with all manual grinders & polishers. The dosing arm with dosing nozzles may be freely positioned on the grinding & polishing machine. The dosing is controlled via the Qdoser GO operating panel. The pre-dosing and backflushing function may also be started manually. After the Qdoser GO has been activated, the liquids are dispensed/dosed automatically during the entire preparation procedure. In combination with other QATM grinders & polishers, the dosing is also started automatically without the need for manual activation.

- | Automatic dosing of
 - 3x diamond suspensions and 1x lubricant or
 - 4x diamond suspension Dia-Complete
- | Dosing intervall adjustable
- | Reverse rinsing function
- | Single or automatic mode
- | For all grinders & polishers
- | Glass bottles for suspensions refillable and removable

DOSING BOTTLE FOR GRINDING & POLISHING MACHINES

QDOSER ONE



The Qdoser ONE dosing system is an easy-to-use solution for manual or semi-automatic polishing processes. QATM's standard suspension and lubricant bottles can be quickly screwed on the dosing attachment. The dosing volume is closely and continuously controlled with an adjusting screw. The system can be mounted to the polishing heads Qpol GO and ECO. If used with manual grinding/polishing machines, a stand with suction cups can be placed on the machines.

- | Dosing attachment for standard and suspension bottles
- | Adjustable dosing volume
- | Can be freely positioned with stand
- | Can be screwed to polishing heads Qpol GO and ECO

