

Qprep GALAXY diamond grinding disc

QPREP GALAXY diamond grinding discs are with specially arranged elements, which contains resin-bonded diamond grains, on a stainless-steel metal carrier equipped. They are used for planar and fine grinding of medium-hard and hard materials. The color-coded grinding discs cover the FEPA grain sizes P80 to P1000. The color coding for the individual grain size ranges can be found in the application table.



PRODUCT ADVANTAGES

- | High stock removal
- | Short processing times
- | High planarity

RECOMMENDED APPLICATIONS

- | Steel with medium and high hardness
- | Cast iron (CJS / CJL)
- | Composites
- | Hard metal
- | Glass



Item No.	Unit	Description	
GALAXY DIAMOND GRINDING DISCS BOND: RESIN			
		Description	Ø
for medium to hard materials			
• planar grinding: grey (P80-100)			
• Fixation system: Magnetic foil			
95005505	1 Pc.	GALAXY diamond grinding disc grey	200 mm
95004279	1 Pc.	GALAXY diamond grinding disc grey	250 mm
95004280	1 Pc.	GALAXY diamond grinding disc grey	300 mm
95013917	1 Pc.	GALAXY diamond grinding disc grey	350 mm
incl. dressing stone No. 95006603			
for medium to hard materials			
• planar grinding: red (P100-120)			
• Fixation system: Magnetic foil			
95005521	1 Pc.	GALAXY diamond grinding disc red	200 mm
95004310	1 Pc.	GALAXY diamond grinding disc red	250 mm
95004314	1 Pc.	GALAXY diamond grinding disc red	300 mm
95013918	1 Pc.	GALAXY diamond grinding disc red	350 mm
incl. dressing stone No. 95006603			
for medium to hard materials			
• fine grinding: green (P180-320)			
• Fixation system: Magnetic foil			
95005522	1 Pc.	GALAXY diamond grinding disc green	200 mm
95004311	1 Pc.	GALAXY diamond grinding disc green	250 mm
95004315	1 Pc.	GALAXY diamond grinding disc green	300 mm
95013919	1 Pc.	GALAXY diamond grinding disc green	350 mm
incl. dressing stone No. 95006603			

Item No.	Unit	Description
----------	------	-------------

GALAXY DIAMOND GRINDING DISCS
BOND: RESIN



Description	∅
-------------	---

for medium to hard materials
 • fine grinding: blue (P400-600)
 • Fixation system: Magnetic foil

95005523	1 Pc.	GALAXY diamond grinding disc blue	200 mm
95004312	1 Pc.	GALAXY diamond grinding disc blue	250 mm
95004316	1 Pc.	GALAXY diamond grinding disc blue	300 mm
95013920	1 Pc.	GALAXY diamond grinding disc blue	350 mm

incl. dressing stone No. 95006603

for medium to hard materials
 • fine grinding: yellow (P800-1000)
 • Fixation system: Magnetic foil

95005524	1 Pc.	GALAXY diamond grinding disc yellow	200 mm
95004313	1 Pc.	GALAXY diamond grinding disc yellow	250 mm
95004317	1 Pc.	GALAXY diamond grinding disc yellow	300 mm
95013921	1 Pc.	GALAXY diamond grinding disc yellow	350 mm

incl. dressing stone No. 95006603

Notes
