

## **Correp** Silicon carbide grinding paper with foil backing

Silicon carbide grinding paper with foil backing can be used individually and are ideally suited for a wide range of varied materials. They offer a high degree of flexibility, especially for frequently changing material variants. Furthermore, the foiled paper enables the use of the "Vakujet" vacuum fixation in the QATM automatic grinding and polishing machine Qpol 300 Bot. QPREP silicon carbide grinding paper with foil backing is available in diameters of 200 / 250 / 300 / and 350 mm and in 12 different grain sizes.



### **PRODUCT ADVANTAGES**

- Suitable for vacuum fixation
- Directionally oriented applied abrasive grains with homogeneous distribution for a high cutting effect and efficient material removal
- I Ideally suited for the QPREP adhesive carrier disc
- I 12 grain sizes for a wide range of applications

### **RECOMMENDED APPLICATIONS**

- I All-purpose usage
- Suitable for planar, pre- and fine grinding steps due to different grain sizes

		WITH FOIL BACKING	unio EEDA atausi III
			rain FEPA standard
		<ul> <li>Fixation system: Adhesive carrier disc</li> </ul>	
		200 mm Ø	
95011898	100 Pcs.	SiC grinding paper with foil backing	P80
95011899	100 Pcs.	SiC grinding paper with foil backing	P120
95011900	100 Pcs.	SiC grinding paper with foil backing	P180
95011901	100 Pcs.	SiC grinding paper with foil backing	P240
95011902	100 Pcs.	SiC grinding paper with foil backing	P320
95011903	100 Pcs.	SiC grinding paper with foil backing	P400
95011904	100 Pcs.	SiC grinding paper with foil backing	P500
95011905	100 Pcs.	SiC grinding paper with foil backing	P600
95011906	100 Pcs.	SiC grinding paper with foil backing	P800
95011907	100 Pcs.	SiC grinding paper with foil backing	P1000
95011908	100 Pcs.	SiC grinding paper with foil backing	P1200
95011909	100 Pcs.	SiC grinding paper with foil backing	P2500
		250 mm Ø	
95011928	100 Pcs.	SiC grinding paper with foil backing	P80
95011929	100 Pcs.	SiC grinding paper with foil backing	P120
95011930	100 Pcs.	SiC grinding paper with foil backing	P180
95011931	100 Pcs.	SiC grinding paper with foil backing	P240
95011932	100 Pcs.	SiC grinding paper with foil backing	P320
95011933	100 Pcs.	SiC grinding paper with foil backing	P400
95011934	100 Pcs.	SiC grinding paper with foil backing	P500
95011935	100 Pcs.	SiC grinding paper with foil backing	P600
95011936	100 Pcs.	SiC grinding paper with foil backing	P800
95011937	100 Pcs.	SiC grinding paper with foil backing	P1000
95011938	100 Pcs.	SiC grinding paper with foil backing	P1200
95011939	100 Pcs.	SiC grinding paper with foil backing	P2500









# SILICON CARBIDE GRINDING PAPER WITH FOIL BACKING

### **Grain FEPA standard**

• Fixation system: Adhesive carrier disc

### 300 mm Ø

95	011943	100 Pcs.	SiC grinding paper with foil backing	P80
95	011944	100 Pcs.	SiC grinding paper with foil backing	P120
95	011945	100 Pcs.	SiC grinding paper with foil backing	P180
95	011946	100 Pcs.	SiC grinding paper with foil backing	P240
95	011947	100 Pcs.	SiC grinding paper with foil backing	P320
95	011948	100 Pcs.	SiC grinding paper with foil backing	P400
95	011949	100 Pcs.	SiC grinding paper with foil backing	P500
95	011950	100 Pcs.	SiC grinding paper with foil backing	P600
95	011951	100 Pcs.	SiC grinding paper with foil backing	P800
95	011952	100 Pcs.	SiC grinding paper with foil backing	P1000
95	011953	100 Pcs.	SiC grinding paper with foil backing	P1200
95	011954	100 Pcs.	SiC grinding paper with foil backing	P2500

		350 mm Ø	
95014005	100 Pcs.	SiC grinding paper with foil backing	P80
95014006	100 Pcs.	SiC grinding paper with foil backing	P120
95014007	100 Pcs.	SiC grinding paper with foil backing	P180
95014008	100 Pcs.	SiC grinding paper with foil backing	P240
95014009	100 Pcs.	SiC grinding paper with foil backing	P320
95014010	100 Pcs.	SiC grinding paper with foil backing	P400
95014011	100 Pcs.	SiC grinding paper with foil backing	P500
95014012	100 Pcs.	SiC grinding paper with foil backing	P600
95014013	100 Pcs.	SiC grinding paper with foil backing	P800
95014014	100 Pcs.	SiC grinding paper with foil backing	P1000
95014015	100 Pcs.	SiC grinding paper with foil backing	P1200
95014016	100 Pcs.	SiC grinding paper with foil backing	P2500

Ν	0	tε	es	