



| Item No. | Unit | Description | |
|---|--------|-------------|--------|
| POLISHING CLOTHS, SELF-ADHESIVE | | | |
| Description Ø | | | |
| ZETA – Polishing cloth | | | |
| <ul style="list-style-type: none"> • Suitable for end polishing as well as fine polishing with fine polishing suspension • Short flocked, soft synthetic cloth • for all materials • For use with diamond grain size: 3/1 µm • Suitable for oxide suspensions as well as alumina • Fixation system: X-Tap | | | |
| 92005681 | 5 Pcs. | Zeta | 200 mm |
| 92005683 | 5 Pcs. | Zeta | 250 mm |
| 92005684 | 5 Pcs. | Zeta | 300 mm |
| 95016037 | 5 Pcs. | Zeta | 350 mm |
| KAPPA – Polishing cloth | | | |
| <ul style="list-style-type: none"> • Suitable for fine polishing with fine polishing suspension • Soft polishing felt • For all materials • Suitable for alumina • Fixation system: X-Tap | | | |
| 92002563 | 5 Pcs. | Kappa | 200 mm |
| 92002572 | 5 Pcs. | Kappa | 250 mm |
| 92002581 | 5 Pcs. | Kappa | 300 mm |
| 95016038 | 5 Pcs. | Kappa | 350 mm |
| OMEGA – Polishing cloth | | | |
| <ul style="list-style-type: none"> • Suitable for fine polishing with fine polishing suspension • Soft synthetic cloth (chemical resistant) • For all materials • Suitable for oxide suspensions • Fixation system: X-Tap | | | |
| 92002562 | 5 Pcs. | Omega | 200 mm |
| 92002571 | 5 Pcs. | Omega | 250 mm |
| 92002580 | 5 Pcs. | Omega | 300 mm |
| 95016039 | 5 Pcs. | Omega | 350 mm |
| LAMBDA – Polishing cloth | | | |
| <ul style="list-style-type: none"> • Suitable for fine polishing with fine polishing suspension • Brushed, structured synthetic cloth (chemical resistant) • For all materials • Suitable for oxide suspensions as well as alumina • Fixation system: X-Tap | | | |
| 95003597 | 5 Pcs. | Lambda | 200 mm |
| 95003598 | 5 Pcs. | Lambda | 250 mm |
| 95003599 | 5 Pcs. | Lambda | 300 mm |
| 95016040 | 5 Pcs. | Lambda | 350 mm |

all polishing cloths also available for Qpol 30 (Ø 30 mm)