

IR100FD Series Automatic Refractometers (Fully Compliant 21CFR ,GLP ,GMP)



Enjoy the latest FD+ software

Specification		Model	IR120FD+	IR140FD+	IR180FD+
Measurement parameters		Same as IR100 series , refer to Page 12			
Calibration		1 point calibration		1 point calibration or up to 5 points calibration	
21CFR Electronic Features	User Management	No. of users	50		
		User entry password	Yes		
		User level	4		
	Data management	Data processing	PW protection		
		PW reminder	Yes		
		Measuring mode	200,editable		
		Discarded method database	1000 Groups data		
		Filed measurement database	100000 Groups data		
	Safe guard in data transmission	Discarded measurement database	10000 Groups data		
		Audit trail database	1000000 Groups data		
Software to PC with protector		Provided			
Additional features for 21CFR	Data format	PDF ,excel,MD5 verify			
	On-line printing	Yes			
		Database printing	Yes		
Ordering Information					
N.W. / G.W.		6kg / 7kg			
Dimension(L*W*H)		36*24*15cm(instrument),25*16*14cm (printer)			
With packing dimension		51*40*29cm(instrument),39*22*22cm (printer)			
Configuration		IR120FD+ ,Printer Data transfer safety SW Power cord from Page 31	IR140FD+ ,Printer Data transfer safety SW Power cord from Page 31	IR180FD+ ,Printer Data transfer safety SW Power cord from Page 31	
Order number		85001001	85001002	85000055	