

Homogenizer Digital Model

NEW!

270



High Safety

Overheat-Heat and Current Protection function prevent malfunction of motor.


Safe & Fine Compatibility in Tools

Useful for All Dispersing and Homogenizing. Applications and realize fine homogenizing. High Fidelity FND Display Easy-to-read backlit LCD for display.

Convenient Speed Set Knob Easy to use Speed Set Knob for quick and convenient selection of mixing speed.

Model	SH-HZD		
Speed Range	2,000 to 27,000rpm(without load)		
Maximum Volume (H2O)	2,500ml		
Motor Rating(Input/Output)	400 / 300 W		
Speed Control	Stepless		
Display	3 Digit FND		
Safety Devices	Over Heat Protector		
Permissible Temp. & Humidity	5°C~40°C, <80%		
Usable Dispersing Tool	Ø8mm Stator "SH-HZ1008", Ø10mm Stator "SH-HZ1010 Ø18mm Stator "SH-HZ1018", Ø25mm Stator "SH-HZ1025		
Material	ABS & Aluminum Stainless Steel(#304)-Tool		
Extension Arm Size(ØxL)	12.7x250mm / 0.5x9.84''		
Power Consumption	400W		
Power Supply 1 Phase, AC 120V, 60Hz or AC 230V, 50/60 HZ			

Accessories



SH-HZ1010

SH-HZ1018

SH-HZ1025

SH-OSSR010	Standard Plate Stand, for Overhead Stirrer and Homogenizer, 25x37cm Plate, with Rod and Silicon Pad
SH-OSHC010	Clamp for Overhead Stirrers and Homogenizers, up to Ø23mm Rod

Homogenizer Tool, Rotor ¢7.5, Stator ¢10, 1~100 ml

Homogenizer Tool, Rotor Φ12.7, Stator Φ18, 10~1000 ml

Homogenizer Tool, Rotor Φ17, Stator Φ25, 50~2500 ml